



Gastorr FG-42

Solvent Degassing Unit Gastorr FG-42

Degassing unit for fluorinated solvents

Features

- Meant for use with fluorinated solvents
- Self-cleaning function protects device interior
- CPU control achieves stable degassing
- Independent chamber design prevents solvent contamination

Product Name	Gastorr Solvent Degassing Unit	
Product Number	FG-42	
Channels	2	
Internal Volume	10 mL	
Flow Rate	≤ 5 mL/min	
Vacuum Range	150 hPa fixed	
Thread Size	1/4-28 UNF	
Wetted Material	PTFE, ETFE, PPS	
Remote Control	Input: — Output: Error	
Dimensions	108 (H) × 70 (W) × 240 (D) mm (excluding protuberances)	
Weight	Approx. 2.3 kg	
Power	AC 100 ~ 240 V (50 / 60 Hz) , MAX 80 VA	
Order Code	SGA-G3008	



Gastorr PG-41P

Gastorr PG-42P

Solvent Degassing for Fractionating Applications Gastorr PG Series

Degassing units with flow rates up to 50 mL/min for fractionating applications

Features

- Flow rates up to a max of 50 mL/min
- Self-cleaning function protects device interior
- CPU control achieves stable degassing
- Independent chamber design prevents solvent contamination
- Choose from 1 channel or 2 channel configurations
- Real time OLED display of chamber vacuum level

Product Name	Gastorr Solvent Degassing Unit	
Product Number	PG-41P	PG-42P
Channels	1	2
Internal Volume	62.0 mL	
Flow Rate	≤ 50 mL/min	
Vacuum Range	150 hPa fixed	
Thread Size	1/4-28 UNF	
Wetted Material	PTFE, ETFE, PEEK	
Remote Control	Input: — Output: Error	
Dimensions	180 (H) × 90 (W) × 400 (D) mm (excluding protuberances)	
Weight	Approx. 5.0 kg	Approx. 6.0 kg
Power	AC 100 ~ 240 V (50 / 60 Hz) , MAX 80 VA	
Order Code	SGA-G3501	SGA-G3502