


Multi-Functional Piston Pump <OEM> Switching Delivery Pump DP-11




Dimensions




Product Specifications	
Part Number	DP-11A
Pump Flow Modes	Measuring volume, Continuous flow, One-shot flow
Numbers of ports	6
Measuring Vol. Settings	0.001 ~ 9999.999 mL
Flow Rate Ranges	0.1 ~ 25.0 mL/min
Communication	RS232C, RS485 (between FLOM devices)
Power	DC 24 V , MAX 30 W

Note: See page 10 for the standalone version of this product.

Solvent Degassing Units <OEM> Gastorr AG / BG Series



Dimensions



OEM UNIT

Product Specifications	
Product Number	AG / BG-4xA
Channels	1 ~ 6 (pictured is 4-channel unit)
Wetted Material	Teflon™ AF *, PFA, PPS
Internal Volume	340 µL or 510 µL
Remote Control	Input: —, Output: Error
Display	AG: —, BG: OLED Display
Power	DC 24 V , MAX 30 W

*:Teflon™ AF is a registered trademark of the Chemours Corporation

Note: See page 17 for the standalone version of this product.

OLED Display

