



7. Other Components

○ Filters

Solvent Filters




Name: Suction Filter
P/N: 8810
Pore Size: 20 μm
Connection: OD 1/8", 3.0 mm Tubing
Material: PTFE, PCTFE
Quantity: 1




Name: Suction Filter
P/N: 8811
Pore Size: 20 μm
Connection: OD 1/16" Tubing
Material: PTFE, PCTFE
Quantity: 1

Solvent Filters (Inline)



Name: PEEK Inline Filter
P/N: 8502
Pore Size: 5 μm
Connection: No.10-32 UNF
Material: PEEK
Quantity: 1



Name: PEEK Inline Filter
P/N: 8501
Pore Size: 10 μm
Connection: No.10-32 UNF
Material: PEEK
Quantity: 1



Name: PEEK Frit
(Inline Filter Replacement)
P/N: 8512
Pore Size: 5 μm
Material: PEEK
Quantity: 5



Name: PEEK Frit
(Inline Filter Replacement)
P/N: 8511
Pore Size: 10 μm
Material: PEEK
Quantity: 5

Note: The dimensions listed are in mm.

○ Solvent Bottles



Name: Solvent Bottle
P/N: 9841
Capacity: 1000 mL
Cap Type: 9851 Bottle Cap
Quantity: 1



Name: Solvent Bottle
P/N: 9842
Capacity: 500 mL
Cap Type: 9851 Bottle Cap
Quantity: 1



Name: Bottle Cap
P/N: 9851
Type: GL45 (for Schott Solvent Bottles)
Tubing OD can be 1/16", 1/8", 3.0 mm, 2.0 mm
Quantity: 1

○ Tube Cutters



Name: PEEK Tube Cutter
P/N: 9830
Quantity: 1



Name: SUS Tube Cutter
P/N: 9831
Quantity: 1