


○ETFE Fittings and Joints for OD 1/16", 1/8", 3.0 mm Tubing

ETFE Fitting and Ferrule (No.10-32 UNF, OD 1/16" Tubing)




Name: 1/16" Fitting D
P/N: 9820D
Material: ETFE
Quantity: 10/pack




Name: 1/16" Ferrule D
P/N: 9821D
Material: ETFE
Quantity: 10/pack


ETFE FLAT SEAL Series (1/4-28 UNF, OD 3.0 mm or 1/8" Tubing)




Name: 1/8" HEX Flat Seal Fitting D
P/N: 9822D
Material: ETFE
Quantity: 10/pack



Name: 1/8" Flat Seal Fitting D
P/N: 9823D
Material: ETFE
Quantity: 10/pack




Name: 1/8" Flat Seal Ferrule
P/N: 9774
Material: ETFE
Quantity: 10/pack




Name: 3.0 mm Flat Seal Ferrule
P/N: 9775
Material: ETFE
Quantity: 10/pack


ETFE T-joints



Name: T-joint D
P/N: 9824D
Material: ETFE
Thread: No.10-32 UNF
ID: 1.0 mm
Quantity: 1



Name: T-joint D
P/N: 9825D
Material: ETFE
Thread: 1/4-28 UNF
ID: 1.0 mm
Quantity: 1



Name: Multi-thread T-joint D
P/N: 9826D
Material: ETFE
Thread: No.10-32 UNF, 1/4-28 UNF
ID: 1.0 mm
Quantity: 1