



○ **Other Polymer Fittings for OD 1/16", 1/8", 3.0 mm, and 2.0 mm Tubing**
Low Pressure Fittings (up to 1 MPa)


FLAT SEAL Series (1/4-28UNF, OD 1/16" or 2.0 mm Tubing)




Name: 2.0 mm HEX Flat Seal Fitting
P/N: 9772
Material: PEEK
Quantity: 10/pack




Name: 2.0 mm M6 Flat Seal Ferrule (Penetrating)
P/N: 9755
Material: ETFE
Quantity: 10/pack



Name: 2.0 mm M6 Flat Seal Ferrule (Non-Penetrating)
P/N: 9705
Material: ETFE
Quantity: 10/pack




Name: 1/16" HEX Flat Seal Fitting
P/N: 9773
Material: PEEK
Quantity: 10/pack




Name: 1/16" M6 Flat Seal Ferrule
P/N: 9706
Material: ETFE
Quantity: 10/pack


FLAT SEAL Series (1/4-28UNF, OD 1/8" or 3.0 mm Tubing)



Name: 1/8" HEX Flat Seal Fitting
P/N: 9771
Material: PEEK
Quantity: 10/pack




Name: 1/8" Flat Seal Ferrule
P/N: 9774
Material: ETFE
Quantity: 10/pack




Name: 3.0 mm Flat Seal Ferrule
P/N: 9775
Material: ETFE
Quantity: 10/pack

FLAT SEAL Plug (1/4-28 UNF)



Name: Flat Seal Plug
P/N: 9720P
Material: PP
Quantity: 10/pack

Stainless Steel Fitting (M8×1.25, OD 1/8" Tubing)




Name: 1/8" Fitting
P/N: 9306
Material: SUS316
Quantity: 10/pack




Name: 1/8" Ferrule
P/N: 9321
Material: SUS316
Quantity: 10/pack


DOUBLE LOCK Series (M8×1.25, OD 3.0 mm or 1/8" Tubing)




Name: 1/8" Double Lock Fitting
P/N: 9303
Material: SUS316
Quantity: 10/pack



Name: 1/8" Double Lock Fitting
P/N: 9303P
Material: PEEK
Quantity: 10/pack



Name: 1/8" Double Lock Ferrule
P/N: 9402
Material: ETFE
Quantity: 10/pack



Name: 3.0 mm Double Lock Ferrule
P/N: 9405
Material: ETFE
Quantity: 10/pack

Note: The dimensions listed are in mm.