



○ Polymer Unions and Joints

Unions




Name: Union
P/N: 9005
Material: PEEK
Thread: No.10-32 UNF
ID: 0.5 mm
Quantity: 1




Name: Union
P/N: 9101
Material: PEEK
Thread: 1/4-28 UNF
ID: 1.0 mm
Quantity: 1


Reducing Unions



Name: Reducing Union
P/N: 9110
Material: PEEK
Thread: 1/4-28 UNF, No.10-32 UNF
ID: 1.0 mm
Quantity: 1




Name: Reducing Union
P/N: 9330P
Material: PEEK
Thread: M8×1.25, No.10-32 UNF
ID: 1.0 mm
Quantity: 1




Name: Reducing Union
P/N: 9112
Material: PEEK
Thread: No.10-32 UNF, Syringe
ID: 1.0 mm
Quantity: 1


T-joints



Name: T-joint
P/N: 9006
Material: PEEK
Thread: No.10-32 UNF
ID: 0.5 mm
Quantity: 1




Name: T-joint
P/N: 9131
Material: PEEK
Thread: 1/4-28 UNF
ID: 1.0 mm
Quantity: 1




Name: Multi-thread T-joint
P/N: 9108
Material: PEEK
Thread: No.10-32 UNF, 1/4-28 UNF
ID: 1.0 mm
Quantity: 1

Cross Joints



Name: Cross joint
P/N: 9007
Material: PEEK
Thread: No.10-32 UNF
ID: 0.5 mm
Quantity: 1



Name: Cross joint
P/N: 9109
Material: PEEK
Thread: 1/4-28 UNF
ID: 1.0 mm
Quantity: 1

Note: The dimensions listed are in mm.