


6. Tubing and Fittings

○Stainless Steel Fittings for 1/16” Tubing
Ultra-high Pressure Fittings (up to 140 MPa)



Name: Dual-Fitt SET
P/N: 9620
Material: SUS316, PEEK
Quantity: 10/pack



Name: Dual-Fitting
P/N: 9621
Material: SUS316
Quantity: 10/pack



Name: Replacement Ferrule
P/N: 9623-05 (5/pack), 9623-10 (10/pack)
Material: SUS316, PEEK
Quantity: 5/pack, 10/pack

Dual-Fitt

Reuseable high pressure joints

Features

- Reduces tube twisting
- Ferrule does not bite into tube, making reuse easy

Tightening guide


Part No.	Pressure	Sealing Torque	Tube
9620	140 MPa (20000 psi)	1.4 N•m	Stainless Steel tube
	100 MPa (14500 psi)	0.8 N•m	
	50 MPa (7250 psi)	0.4 N•m	

*The guide is just for reference purposes.

Medium Pressure Fittings (up to 40 MPa)
SUS Ferrule Type (No.10-32 UNF, OD 1/16” Tubing)



Name: 1/16” Fitting S
P/N: 9501
Material: SUS316
Quantity: 10/pack



Name: 1/16” Fitting L
P/N: 9502
Material: SUS316
Quantity: 10/pack



Name: 1/16” Ferrule
P/N: 9601
Material: SUS316
Quantity: 10/pack

Note: The dimensions listed are in mm.