


○PEEK Fittings for 1/16" Tubing

Medium Pressure Fittings (up to PEEK: 30 MPa)



Single Lock Fitting for 1/16" Tubing


An unique fitting which prevents tube twisting and ferrule falling down while tightening




Tightening guide

Pressure	Sealing Torque	Tube
30 MPa (4350 psi)	0.3 N•m	PEEK tube

*The guide is just for reference purposes.



Name: 1/16" HEX Single Lock Fitting P
P/N: 9903P
Material: PEEK
Quantity: 10/pack




Name: 1/16" HEX Single Lock Fitting Mini P
P/N: 9935P
Material: PEEK
Quantity: 10/pack


Medium Pressure Fittings (20 ~ 40 MPa)

EASY FITT Series (No.10-32 UNF, OD 1/16" Tubing)


Sealing Torque: 0.3 N•m (at pressure of 40 MPa)
Note: when using OD 1/16" ID 0.13 mm PEEK Tubing




Name: Easy Fitt 10 mm
P/N: 9001
Material: PEEK
Quantity: 10/Pack



Name: Easy Fitt 12 mm
P/N: 9002
Material: PEEK
Quantity: 10/Pack




Name: Easy Fitt Type 2
P/N: 9130
Material: PEEK
Quantity: 10/Pack




Name: Stop Fitt
P/N: 9004
Material: PEEK
Quantity: 10/Pack

HEX FITT Series (No.10-32 UNF, OD 1/16" Tubing)



Name: Hex Fitt
P/N: 9003
Material: PEEK
Quantity: 10/Pack



Name: Hex Fitt Mini
P/N: 9135
Material: PEEK
Quantity: 10/Pack

Note: The dimensions listed are in mm.