



○Stainless Steel Unions, Joints and Tubing


Stainless Steel Unions



Name: Union
P/N: 9511
Material: SUS316
Thread: No.10-32 UNF
ID: 0.5 mm
Quantity: 1




Name: 1/8"×1/16" Reducing Union
P/N: 9309
Material: SUS316
Thread: M8×1.25, No.10-32 UNF
ID: 1.0 mm
Quantity: 1




Name: 1/8" Union
P/N: 9330
Material: SUS316
Thread: M8×1.25
ID: 2.0 mm
Quantity: 1


Stainless Steel T-joints (for OD 1/16" Tubing)



Name: T-joint 0.3
P/N: 9513-03
Material: SUS316
Thread: No.10-32 UNF
ID: 0.3 mm
Quantity: 1




Name: T-joint 0.5
P/N: 9513-05
Material: SUS316
Thread: No.10-32 UNF
ID: 0.5 mm
Quantity: 1




Name: T-joint 1.0
P/N: 9513-10
Material: SUS316
Thread: No.10-32 UNF
ID: 1.0 mm
Quantity: 1

Stainless Steel T-joint (for OD 3.0 mm or 1/8" Tubing)




NEW
Name: 1/8" T-joint
P/N: 9314
Material: SUS316
Thread: M8×1.25
ID: 2.0 mm
Quantity: 1

Stainless Steel Cross Joint (for OD 1/16" Tubing)




Name: Cross Joint 1.0
P/N: 9520
Material: SUS316
Thread: No.10-32 UNF
ID: 1.0 mm
Quantity: 1


Stainless Steel Tubing (OD 1/16" or 1/8")




Name: SS Tubing
P/N: 7021
Material: SUS316
OD: 1/16"
ID: 0.25 mm
Length: 10 m




Name: SS Tubing
P/N: 7022
Material: SUS316
OD: 1/16"
ID: 0.5 mm
Length: 10 m



Name: SS Tubing
P/N: 7023
Material: SUS316
OD: 1/16"
ID: 0.8 mm
Length: 10 m



Name: SS Tubing
P/N: 7024
Material: SUS316
OD: 1/16"
ID: 1.0 mm
Length: 10 m



Name: SS Tubing
P/N: 7025
Material: SUS316
OD: 1/8"
ID: 2.1 mm
Length: 2 m

Note: The dimensions listed are in mm.