

○ PEEK, PTFE, ETFE, PFA Tubing

PEEK Tubing (OD 1/16")



Name: PEEK Tubing
P/N: 5999
Line Color: Red
OD: 1/16"
ID: 0.13 mm
Length: 5 m (P/N: 5999-05)
10 m (P/N: 5999-10)



Name: PEEK Tubing
P/N: 6001
Line Color: Blue
OD: 1/16"
ID: 0.25 mm
Length: 5 m (P/N: 6001-05)
10 m (P/N: 6001-10)



Name: PEEK Tubing
P/N: 6002
Line Color: Orange
OD: 1/16"
ID: 0.50 mm
Length: 5 m (P/N: 6002-05)
10 m (P/N: 6002-10)




Name: PEEK Tubing
P/N: 6003
Line Color: Green
OD: 1/16"
ID: 0.75 mm
Length: 5 m (P/N: 6003-05)
10 m (P/N: 6003-10)




Name: PEEK Tubing
P/N: 60031
Line Color: Gray
OD: 1/16"
ID: 1.0 mm
Length: 5 m (P/N: 60031-05)
10 m (P/N: 60031-10)

PEEK Tubing (OD 1/8")




Name: PEEK Tubing
P/N: 6004A
OD: 1/8"
ID: 1/16"
Length: 10 m




Name: PEEK Tubing
P/N: 60042A
OD: 1/8"
ID: 2.1 mm
Length: 10 m

PTFE Tubing



Name: PTFE Tubing
P/N: 7213
OD: 1/8"
ID: 1/16"
Length: 10 m



Name: PTFE Tubing
P/N: 7005
OD: 3.0 mm
ID: 2.0 mm
Length: 10 m



Name: PTFE Tubing
P/N: 7211
OD: 1/16"
ID: 0.25 mm
Length: 10 m



Name: PTFE Tubing
P/N: 7210
OD: 1/16"
ID: 0.50 mm
Length: 10 m



Name: PTFE Tubing
P/N: 7212
OD: 1/16"
ID: 0.80 mm
Length: 10 m



Name: PTFE Tubing
P/N: 7215
OD: 1/16"
ID: 1.0 mm
Length: 10 m

ETFE Tubing



Name: ETFE Tubing
P/N: 7009
OD: 1/16"
ID: 0.20 mm
Length: 10 m



Name: ETFE Tubing
P/N: 7010
OD: 1/16"
ID: 0.30 mm
Length: 10 m



Name: ETFE Tubing
P/N: 7011
OD: 1/16"
ID: 0.50 mm
Length: 10 m



Name: ETFE Tubing
P/N: 7012
OD: 1/16"
ID: 0.80 mm
Length: 10 m




Name: ETFE Tubing
P/N: 7013
OD: 1/16"
ID: 1.0 mm
Length: 10 m


PFA Tubing



Name: PFA Tubing
P/N: 7103
OD: 2.0 mm
ID: 1.0 mm
Length: 10 m



Name: PFA Tubing
P/N: 7104
OD: 3.0 mm
ID: 2.0 mm
Length: 10 m



Name: PFA Tubing
P/N: 7105
OD: 1/8"
ID: 1/16"
Length: 10 m