

## 4. Liquid Level Sensor



Sonic Reservoir Sensor SR-13

| Product Name    | Sonic Reservoir Sensor                                  |
|-----------------|---------------------------------------------------------|
| Product Number  | SR-13                                                   |
| Wetted Material | Solvent probe: non-contact, Waste probe: PTFE           |
| Remote Control  | Input : —<br>Output: Error                              |
| Communication   | RS232C (option)                                         |
| Dimensions      | 115 (H) × 85 (W) × 80 (D) mm, (excluding protuberances) |
| Weight          | Aprox. 0.7 kg                                           |
| Power           | DC 24V , MAX 15W (AC Adapter standard)                  |
| Order Code      | FA007-500                                               |

### Sonic Reservoir Sensor SR-13

Non-contact liquid level sensor using sound

#### Features

- ☐ Monitors liquid levels in solvent bottles and waste tanks
- ☐ Liquid level displays
- ☐ Alarm sounds at danger level
- ☐ One-touch connection
- ☐ Control with RS232C protocol (option)



#### Probes and Optional Accessories

##### SR-13 Accessories

| P/N       | Name            |
|-----------|-----------------|
| FD007-502 | Solvent Probe   |
| FD007-503 | Waste Probe     |
| FD007-504 | Solvent Bottle  |
| FD007-505 | Waste Tank      |
| FD007-506 | RS232C Terminal |



Solvent Probe



Waste Probe



Solvent Bottle  
(1000 mL, Glass)



Waste Tank  
(20L, PE)