



Compact Manual Valve
V200



Compact Manual Valve
V239

Compact Manual Valves

A convenient, ultra-compact valve for switching flow paths

Features

- ☐ Attach to the valve panel (sold separately) to allow operation from valve head side
- ☐ A line-up of 5 units to cover a wide variety of applications

Product Name	Compact Manual Valve				
Product Number	V200	V201	V204	V512	V239
Type	6 port 2 position	3 port 2 way	2 port selection	2 port selection	2 port selection
Path Position					
Max Pressure	20 MPa (2900 psi)	20 MPa (2900 psi)	20 MPa (2900 psi)	30 MPa (4350 psi)	1 MPa (145 psi)
Wetted Material	PEEK	PEEK	PEEK	SUS316 , PEEK	PEEK
Tube OD	1/16"	1/16"	1/16"	1/16"	3.0 mm or 1/8"
Path Diameter	0.6 mm (0.024")	0.6 mm (0.024")	0.8 mm (0.032")	0.8 mm (0.032")	1.0 mm (0.039")
Thread Size	No.10-32 UNF	No.10-32 UNF	No.10-32 UNF	No.10-32 UNF	1/4-28 UNF
Order Code	VH-V200	VH-V201	VH-V204	VH-V512	VH-V239



Manual Valve V05

Manual Valve

A convenient manual valve for switching between flow paths.

Features

- ☐ Low torque design for easy switching
- ☐ Maximum pressure 35MPa (5075 psi)

Product Name	Manual Valve
Product Number	V05
Type	6 port 2 position
Path Position	
Max Pressure	35 MPa (5075 psi)
Wetted Material	SUS316 , PEEK
Tube OD	1/16"
Path Diameter	0.6 mm (0.024")
Thread Size	No.10-32 UNF
Order Code	VH-V05



Injector Panel

P/N	Name
4515	Injector Panel

* The valve shown in the picture is not included and must be purchased separately.

*This cannot be attached to FLOM pumps. See page 11 for side panels for valve attachment to FLOM pumps.



Compact Valve Panel (V200, 201, 204, 239, V512)

P/N	Name
45150	Valve Panel

* The valve shown in the picture is not included and must be purchased separately.

*This cannot be attached to FLOM pumps. See page 11 for side panels for valve attachment to FLOM pumps.